
























Laboratory Ergonomics	Features/Options	Image	Applications
<p><b>Arm Supports/Padding</b></p> <p>Wedge-Ease Ergonomic Forearm Supports  3" wedge 79597  2" Wedge 79598  Regular Stacker 79599  Mini Stacker 79600</p> <p>Source:  <a href="http://www.alimed.com/ProductDetail.asp?style=79597&amp;fprd=Wedge-Ease+Ergonomic+Forearm+Supports&amp;oid1=&amp;oid2=">http://www.alimed.com/ProductDetail.asp?style=79597&amp;fprd=Wedge-Ease+Ergonomic+Forearm+Supports&amp;oid1=&amp;oid2=</a></p>	<ul style="list-style-type: none"> <li>- Vinyl w/ nonslip base</li> <li>- Four different sizes can be combined/stacked</li> <li>- Also available in ESD</li> <li>- Rated Fed Std. 209 and ISO14644 Clean Room Class 100 III</li> </ul>		<ul style="list-style-type: none"> <li>- Can be used to provide angled and flat forearm support and to protect against contact stress. Can be used in front of microscopes and to support forearms when doing precision tasks. Four styles available in sloped wedge or flat. Can be stacked.</li> </ul>
<p>Morency Forearm support for microscopes</p> <p>Source:  <a href="http://www.morencyrest.com/msm.htm">http://www.morencyrest.com/msm.htm</a>  MSM-N-BL – 8" wide  MSM-N-BL-T – 5" wide</p>	<ul style="list-style-type: none"> <li>- 8" wide and 5" wide versions.</li> <li>- include a clamp the fits from 5/8 to 1 3/4 inch thick work surfaces. Adaptor available for thicker workstations</li> </ul>		<p>Angle-adjustable forearm support for microscope users</p>
<p>Deluxe Edge Rest 73075</p> <p>Source:  <a href="http://www.alimed.com/ProductDetail.asp?style=73075">http://www.alimed.com/ProductDetail.asp?style=73075</a>  Ergomart  Ergoworks</p>	<ul style="list-style-type: none"> <li>- Attaches securely with self-adhesive hook and loop tabs.</li> <li>- 1 1/8" thick expanded urethane</li> <li>- Black. 22"L x 5"W</li> <li>- Self-skinned foam surface wipes clean.</li> </ul>		<ul style="list-style-type: none"> <li>- Fits over top edge of table surface for protection from sharp edges</li> </ul>
<p>Alimed Ulnar Gel Pads 78336</p> <p>Source:  <a href="http://www.alimed.com/ProductDetail.asp?style=78336&amp;fprd=AliMed%26reg%3B+Ulnar+Gel+Pads+and+Hand+Rest+Gel+Sleeves&amp;oid1=&amp;oid2=">http://www.alimed.com/ProductDetail.asp?style=78336&amp;fprd=AliMed%26reg%3B+Ulnar+Gel+Pads+and+Hand+Rest+Gel+Sleeves&amp;oid1=&amp;oid2=</a></p>	<ul style="list-style-type: none"> <li>- Gel pads to place under elbow</li> <li>- Sold in pairs</li> </ul>		<ul style="list-style-type: none"> <li>- Useful for pipetting and other high-precision tasks where the elbow rests on a hard surface</li> </ul>





Laboratory Ergonomics	Features/Options	Image	Applications
<p>AliEdge Edge Protector, 3'L -75841</p> <p>AliEdge Edge Protector, 6'L - 75842</p> <p>Source:  <a href="http://www.alimed.com/ProductDetail.asp?style=75841&amp;fprd=AliEdge%26trade%3B&amp;oid1=&amp;oid2=">http://www.alimed.com/ProductDetail.asp?style=75841&amp;fprd=AliEdge%26trade%3B&amp;oid1=&amp;oid2=</a></p>	<ul style="list-style-type: none"> <li>- Expanded polyethylene edge rest. A hook and loop fastener is included for easy attachment. Closed-cell foam wipes clean 36" long or 72" long x 3 1/2" wide, 1/2" thick. Overhang is 1 3/8"</li> <li>- Also available for corners</li> </ul>		<ul style="list-style-type: none"> <li>- Edge rest.</li> </ul> <p>Can be used on the front edge of biosafety cabinets away from the downdraft</p>
<p>SoftEdge Corners 75189/GREY/NA</p> <p>Source:  <a href="http://www.alimed.com/ProductDetail.asp?style=75189">http://www.alimed.com/ProductDetail.asp?style=75189</a></p>	<p>4/pk</p>		<p>Cover sharp corners. Can be used with Aliedge edge protectors</p>
<p>Custom padded sleeves worn under lab coat</p> <p>Source:  <a href="mailto:lindalorrainegloves@feri.com">lindalorrainegloves@feri.com</a></p> <p>Contact Linda directly to get quote and make arrangements-  T: 415-239-2397  F: 415-239-6231</p>	<ul style="list-style-type: none"> <li>- Custom-made to fit by Linda Lorraine</li> <li>- Light and comfortable</li> <li>- Mesh fabric breathes well.</li> <li>- Need to separately obtain closed-cell foam for pads at elbow and/or forearm (see below)</li> </ul>	 	<p>Wearable pads can provide elbow and forearm protection from hard surfaces/edges.</p>
<p>Closed-cell foam to be used in sleeves above (cut to fit)</p> <p>Source:  <a href="http://www.softknees.com/softknees.html">http://www.softknees.com/softknees.html</a></p>			<p>SoftKnees pads cut for elbows &amp; forearm pads to fit inside sleeves above</p>




Laboratory Ergonomics	Features/Options	Image	Applications
<p>Nuair Foam Elbow Pads</p> <p>Source:  <a href="http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=30">http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=30</a></p>	<ul style="list-style-type: none"> <li>- Disposable white foam pad</li> <li>- 3/16" closed cell foam</li> <li>- Does not absorb liquid</li> <li>- Package of 10</li> <li>- Self-stick adhesive tape</li> </ul>		<p>Easy to remove and replace foam elbow pads used to pad surface of biological safety cabinet to reduce soft tissue compression and discomfort.</p>
<b>Micromanipulation – De-Capping Devices</b>			
<p>Decapitator 2</p> <p>Source:  Drummond Scientific Company</p> <p><a href="http://www.pipetaid.com/pages/new/new04.html">http://www.pipetaid.com/pages/new/new04.html</a></p> <p>Contact: Chuck Locke  800-523-7480</p>	<ul style="list-style-type: none"> <li>- Accommodates tube sizes from 0.2 mL to 2.2 mL</li> <li>- Comes with a sturdy lab bench stand.</li> </ul>		<p>One-handed method for uncapping micro centrifuge tubes. The unit provides two slots into which the user can insert the appropriate sized tube cap. Once secured, a slight downward pressure opens the tube.</p>
<p>MicroTube Tool (Pack of 4) EF4321</p> <p>Source:  <a href="http://www.daigger.com/catalog/product?deptId=Microcentrifuge+Tube+Tools&amp;prodId=4321">http://www.daigger.com/catalog/product?deptId=Microcentrifuge+Tube+Tools&amp;prodId=4321</a></p>	<p>-</p>		<p>MicroTube Tool  Saves wear and tear on your fingers and thumbs!</p> <p>Innovative MicroTube Tool makes it a breeze to open snap caps and screw caps. Convenient two sided, polypropylene tool features snap-cap opener on one side and a screw-cap device on the other. Pack of 4.</p>

Laboratory Ergonomics	Features/Options	Image	Applications
<p>Microcentrifuge Tube Opener 199250001</p> <p>Source: <a href="http://www.belart.com/shop/199250001-microcentrifuge-tube-opener-p-199250001.html?cpath=18_114&amp;oscsid=xogo-eikoq">http://www.belart.com/shop/199250001-microcentrifuge-tube-opener-p-199250001.html?cpath=18_114&amp;oscsid=xogo-eikoq</a></p>	<ul style="list-style-type: none"> <li>- 79mm (3 1/8") long</li> <li>- 3 per bag</li> <li>- One end opens 1.5 to 2.2ml tubes</li> <li>- Forked blade on opposite end opens smaller sizes: 250m, 400m, and 550m tubes</li> </ul>		<p>Tube opener easily open all sizes of micro centrifuge tubes with just one hand.</p>
<b>Jar Openers</b>			
<p>Dycem Multi-Purpose Jar Openers</p> <p><a href="http://www.amazon.com/dp/B002BUKCV4/ref=asc_df_B002BUKCV4845487?smid=A2VSTQ09J8RDTA&amp;tag=dealt15484-20&amp;linkCode=asn">http://www.amazon.com/dp/B002BUKCV4/ref=asc_df_B002BUKCV4845487?smid=A2VSTQ09J8RDTA&amp;tag=dealt15484-20&amp;linkCode=asn</a></p>	<ul style="list-style-type: none"> <li>- Non-slip plastic dome or sheet</li> <li>- Sheet diameter is 5 1/2"</li> <li>- Dome measures 2" high and 4 1/2" in diameter</li> <li>- 1/8" thick</li> <li>- Latex free.</li> </ul>		<p>Facilitate removal of stubborn jar and bottle lids. Simply place nonslip plastic dome or sheet on lid, grasp tightly and twist.</p>
<p>Un-Skru Jar Opener</p> <p>Source: <a href="http://www.dynamic-living.com/product/unskru-jar-opener">http://www.dynamic-living.com/product/unskru-jar-opener</a></p>	<ul style="list-style-type: none"> <li>- Must be attached to underside of cabinet base</li> <li>- Can open jars using one hand</li> <li>-</li> </ul>		<p>Wedge the jar lid firmly between the metal gripper post and the appropriate edge of the jar opener and twist. (The metal post can grip lids well because it has 'teeth'.) The white Un-Skru jar and bottle opener will open bottle caps and lids up to 4.75" in diameter.</p>





Laboratory Ergonomics	Features/Options	Image	Applications
<p>Progressive International Multi-Opener</p> <p>Source:  <a href="http://www.amazon.com/Progressive-GT-2950-International-Multi-Opener/dp/B0008D7GO2/ref=pd_sim_k_1">http://www.amazon.com/Progressive-GT-2950-International-Multi-Opener/dp/B0008D7GO2/ref=pd_sim_k_1</a></p>	<ul style="list-style-type: none"> <li>- Opens 6 different types of seals and lids</li> <li>- Durable and comfortable grip</li> <li>- It has a hidden blade that zips open bags</li> </ul>		<p>Jar and seal opener that can be used with a variety of seals, lids and bags.</p>
<p>Micromanipulation Tools and Padding</p>			
<p>Thermo-Pellets 1 # NC32511 for custom handles</p> <p>Source:  <a href="http://www.ncmedical.com/item_961.html">http://www.ncmedical.com/item_961.html</a></p>	<ul style="list-style-type: none"> <li>- 1 and 3 pound packages</li> </ul>		<p>Low-temperature thermoplastic pellets mold together quickly and easily to adapt handles, splints and equipment. The pellets can be used to build up handles. When heated to 140 degrees F (60 degree C), pellets soften quickly and can be formed into a mass that can be molded and remolded as needed. Heat in water, on a heating plate or in a microwave oven. The thermoplastic will stick to most porous surfaces.</p>
<p>Precision hand tools with foam/molded grips</p> <p>Source:  <a href="http://excelta.com/marketplace/browse.asp?category=396">http://excelta.com/marketplace/browse.asp?category=396</a>   <a href="http://www.timemotion.com/Products/Tweezer-5-SA-ET--Excelta-Tweezers__WAS-03-34021.aspx">http://www.timemotion.com/Products/Tweezer-5-SA-ET--Excelta-Tweezers__WAS-03-34021.aspx</a></p>	<ul style="list-style-type: none"> <li>- Multiple styles of precision hand tools</li> <li>- Foam grips</li> </ul>		<p>Molded and foam grips for precision hand tools manufactured by Excelta. Other tool modifications available with bulk orders.</p>





Laboratory Ergonomics	Features/Options	Image	Applications
<p>Low Force Forceps</p> <p>Source:  <a href="http://www.surgipath.com/us-en/Products/5-1-2-ergonomic-forceps">http://www.surgipath.com/us-en/Products/5-1-2-ergonomic-forceps</a></p>	<ul style="list-style-type: none"> <li>- 5 ½" low force forceps</li> <li>- 7" low force forceps</li> </ul>		<p>Reduces force required to hold specimens in forceps while processing.</p>
<p>Dumont Medical Tweezers with Clamping Ring</p> <p>Source:  <a href="http://www.emsdiasum.com/microscopy/products/tweezers/dumont.aspx">http://www.emsdiasum.com/microscopy/products/tweezers/dumont.aspx</a></p>	<ul style="list-style-type: none"> <li>- Variety of sizes</li> <li>- Clamping ring to lock in closed position</li> </ul>		<p>The clamping ring eliminates the need to apply constant finger pressure for gripping.</p>
<p>Negative-Action Style Tweezers (Self-Closing)</p> <p>Source:  <a href="http://www.emsdiasum.com/microscopy/products/tweezers/dumont.aspx">http://www.emsdiasum.com/microscopy/products/tweezers/dumont.aspx</a></p>	<ul style="list-style-type: none"> <li>- Variety of sizes</li> <li>- Default position is closed rather than open</li> </ul>		<p>Negative-Action Style Tweezers (Self-Closing) eliminate the need to apply constant finger pressure for gripping.</p>
Storage Bins, Racks and Holders			
<p>Double Bin Storage SV-7014</p> <p>Source:  <a href="http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=37">http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=37</a></p>	<ul style="list-style-type: none"> <li>- Flip-top Lids</li> <li>- Slide Resistant Rubber Feet</li> <li>- Each Bin: 6" x 5" x 6"</li> </ul>		<p>For Storing Microtubes, Pipet Tips or any Bulk Items</p> <p>Angled openings</p>
<p>Acrylic Dispensing Bins</p> <p>ML9541  ML9542  ML9540</p> <p>Source:  <a href="http://www.marketlabinc.com/products/details/164">http://www.marketlabinc.com/products/details/164</a></p>	<ul style="list-style-type: none"> <li>- ML9541 Small Bin  - 5.25"W x 9"D x 8.5"H</li> <li>- ML9542 Medium Bin  - 6" x 10"D x 10.5"H</li> <li>- ML9540 Large Bin  - 7"W x 12"D x 10.5"H</li> </ul>		<p>The flip top lids on the top and on the front of these general use dispensing bins keeps supplies clean and dust-free. The lid is designed to have a slight overhang, so you can lift it and pick up the items inside with one hand. Holds tube caps, cotton balls, pipette tips, analyzer cups.</p>




Laboratory Ergonomics	Features/Options	Image	Applications
<p>Horizontal Pipette Holder SV-7020</p> <p>Source:  <a href="http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=74">http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=74</a></p>	<ul style="list-style-type: none"> <li>- Size: 20" W x 5" H x 10" D</li> <li>- Accessible Storage of 1ml - 25ml Pipettes</li> <li>- Easy-to-Clean White Plastic</li> <li>- Lightweight and Durable</li> </ul>		<p>Horizontally Mounted Under Shelf or Cabinet</p>
<p>Vertical Pipette Holder SV-7010</p> <p>Source:  <a href="http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=35">http://www.scientificvisions.com/catalog/product_info.php?cPath=32&amp;products_id=35</a></p>	<ul style="list-style-type: none"> <li>- Accessible Storage of 1 ml to 25 ml Pipettes</li> <li>- Easy-to-clean White Plastic</li> <li>- Lightweight and Durable</li> </ul>		<p>Can be Vertically Mounted on Side of Safety Cabinet or Laminar Airflow Hood</p>
<p>Pipette Racks</p> <p><a href="http://www.marketlabinc.com/products/details/1036">http://www.marketlabinc.com/products/details/1036</a></p>	<ul style="list-style-type: none"> <li>- Measures 3.5' w x 11.5'd x 16"h.</li> </ul>	 <p>ML2170</p>	<p>Tilted acrylic pipette rack with four angled compartments for easy retrieval of pipettes. Constructed of durable acrylic. Holds 0.1 to 25ml pipettes or cans and has small footprint.</p>
<p>Turntable Workstations</p> <p>EF1120A – 1 tier  EF1120B – 2 tiers  EF1120 – 3 tiers</p> <p>Source:  <a href="http://www.daigger.com/catalog/product?deptId=&amp;prodId=1120#">http://www.daigger.com/catalog/product?deptId=&amp;prodId=1120#</a></p>	<ul style="list-style-type: none"> <li>- Solid, chemical-resistant PVC, easy 360° rotation. Available with one, two or three tiers. Each tier has four 6 x 6" compartments.</li> </ul>	 <p>2 tiers</p>	<p>Single or stacked turntable organizer can be used in deep cabinets, corners, or shared workstations to reduce reach and optimize storage</p>





Laboratory Ergonomics	Features/Options	Image	Applications
		 <p style="text-align: center;">3 tiers</p> <p style="text-align: center;">1 tier</p>	
<p><b>ML9970 Deluxe Rotating Organizers</b></p> <p>Source:  <a href="http://www.marketlabinc.com/products/details/332">http://www.marketlabinc.com/products/details/332</a></p>	<ul style="list-style-type: none"> <li>- 14"L x 13.5"W x 12.5"H</li> <li>- Easy fingertip rotation for quick access to supplies</li> <li>- Extra large capacity within a small 1 square foot footprint</li> <li>- Built-in dispensers and compartments fit a wide variety of office and lab supplies</li> <li>- Constructed of durable PVC</li> </ul>	 <p style="text-align: right;">ML9970</p>	<p>Turntable organizer can be used in corners or shared workstations to reduce reach and optimize storage</p>
<p><b>Organization Station ML7324</b></p> <p>Source:  <a href="http://www.marketlabinc.com/products/details/469">http://www.marketlabinc.com/products/details/469</a></p>	<ul style="list-style-type: none"> <li>- 14"W x 11"D x 11"H</li> <li>- Customizable microscope organizer</li> <li>- Detached divider boxes allow you to customize your storage; simply attach them with double-sided foam tape included</li> <li>- Constructed of acrylic</li> </ul>	<p><b>ML7324</b></p> 	<p>Customizable microscope organizer with detachable boxes that can be placed where desired with double sided foam.</p>











Laboratory Ergonomics	Features/Options	Image	Applications
<p>ML2117 Medium Step Shelf</p> <p>Source:  <a href="http://www.marketlabinc.com/products/details/344">http://www.marketlabinc.com/products/details/344</a></p>	<ul style="list-style-type: none"> <li>- 18"L x 6"W x 10"H</li> <li>- Utilizes vertical space</li> <li>- Designed to fit over gas and water outlets or other protrusions</li> <li>- Cords and wire fit through space at lower back</li> <li>-</li> </ul>		<p>Instant Step Shelves make use of vertical space, and keep items up and out of the way on wet surfaces. Designed to fit over gas and water outlets or other protrusions. Cords and wires fit through the space at the lower back.</p>
<p>Cleanware, Bag, HDPE, Laboratory Waste, 8X10", 50/Pkg 60-131741008</p> <p>Cleanware, Bag, HDPE, Laboratory Waste, 12x16", 50/Pkg 60-131741012</p> <p>Source:  <a href="http://www.sciencelab.com/page/S/PVAR/23677/60-131741008">http://www.sciencelab.com/page/S/PVAR/23677/60-131741008</a></p>	<ul style="list-style-type: none"> <li>- White bags are made of High Molecular Weight High Density (HMHD) Polyethylene which is steam autoclavable at 121°C (250°F). Maximum weight of contents should not exceed 1.36 kg (3 lb.).</li> </ul>		<p>Disposable waste bags with an adhesive strip across the top frees bench or hood space when hung from walls, cabinets, or virtually anywhere. Can reduce reaching. Peel off the backing and hang close to your work area. The back panel is longer than the front for easy opening.</p>
<p>Turntable NU-977-001</p> <p>Source:  <a href="http://www.scientificvisions.com/catalog/default.php?cPath=26&amp;sort=2a&amp;page=3">http://www.scientificvisions.com/catalog/default.php?cPath=26&amp;sort=2a&amp;page=3</a></p>	<ul style="list-style-type: none"> <li>- 100% Stainless Steel</li> <li>- Autoclavable</li> <li>- Easily Cleaned and Sterilized</li> <li>- 12" Diameter</li> <li>- Stainless Steel Ball Bearing Design</li> <li>-</li> </ul>		<p>Stainless steel rotating table brings work closer to the user to decrease reaching and allows more efficient use of space. Can be used in biosafety cabinets.</p>
<b>Microscope Positioners/Adapters</b>			
<p>Positioning plates \$99.00</p> <p>Source:  <a href="http://www.marketlabinc.com/products/details/11">http://www.marketlabinc.com/products/details/11</a></p>	<ul style="list-style-type: none"> <li>- Solid construction with non-skid surface</li> <li>- Weighted for stability</li> <li>- Positioning Plate = 11.5"L x 8.5"W x 2"H</li> <li>- Height adjustability between 2" to 2.75")</li> <li>-</li> </ul>	 <p style="text-align: center;">ML8350</p>	<p>Good posture is critical when working at a microscope. The Microscope Positioning Plate allows you to adjust your scope's height from 2" to 2.75" by simply dialing the thumb wheels in the corners. The threaded leveling feet on the plate allow you to dial to an infinite number of heights and ergonomic positions.</p>





Laboratory Ergonomics	Features/Options	Image	Applications
<p>Point-O-Care Adjustable Microscope Base AF340</p> <p>Source:  <a href="http://www.medexsupply.com/products/pid-11810/PointOCareAdjustableMicroscope.htm">http://www.medexsupply.com/products/pid-11810/PointOCareAdjustableMicroscope.htm</a></p>	<ul style="list-style-type: none"> <li>- Large 14.5" x 8.5" powder coated steel plate will accept most all brands of microscope</li> <li>- 2" of adjustment from 1" to 3"</li> <li>- Infinite number of heights and angles by adjusting threaded leveling feet</li> <li>- Non-skid feet</li> </ul>		<p>Perfect posture is just a knobs turn away with this adjustable microscope base. Raise and/or tilt your microscope into the ergonomic position you need.</p>
<p>Microscope Stands</p> <p>Source:</p> <ul style="list-style-type: none"> <li>- <a href="http://www.microscopeworld.com/MSWorld/SearchResult.aspx?CategoryID=271&amp;gclid=COWguJmiwJsCFRFMagodmhgHA">http://www.microscopeworld.com/MSWorld/SearchResult.aspx?CategoryID=271&amp;gclid=COWguJmiwJsCFRFMagodmhgHA</a></li> <li>- <a href="http://www.professionalmicroscopes.com/microscope-stands.html">http://www.professionalmicroscopes.com/microscope-stands.html</a></li> </ul>	<ul style="list-style-type: none"> <li>- Available for variety of models and applications</li> <li>- Adjustable height, tilt and horizontal positioning</li> </ul>		<p>Mounting brackets and stands used to position microscopes and microscope components to reduce eye and neck strain associated with prolonged microscope use.</p>
<p>Scopease Microscope Support</p> <p>Source:</p> <ul style="list-style-type: none"> <li>- <a href="http://www.imebinc.com/Item/SC101.htm">www.imebinc.com/Item/SC101.htm</a></li> </ul>	<ul style="list-style-type: none"> <li>- System includes base unit, 2 armrests, 3 tilting bars</li> <li>- Padded armrests</li> <li>- Tilts microscope</li> <li>- 3 models available</li> </ul>		<p>The Scopease system was designed to adjust the position of a variety of microscopes. The soft armrests are adjustable by removing 1" segments of the foam. The microscope can be tilted at any of 9 angles with the adjusting bars.</p>
<p>ErgoAdapter</p> <p>Source:  <a href="http://www.bayoptical.com/ergoadapter.html">http://www.bayoptical.com/ergoadapter.html</a></p>	<ul style="list-style-type: none"> <li>- Simple locking knob adjusts microscope height</li> <li>- High-precision, stainless steel, linear ball bearing</li> <li>- Height adjustment range of +3/4" to +5 inches</li> <li>- Can rotate 180° for left-or right-handed operation</li> <li>- Optical and mechanical alignments constant at any position or inclination angle</li> </ul>		<p>ErgoAdapter inserts between microscope binocular head and optics carrier. When work surfaces are set at elbow height to facilitate hand operations, the ErgoAdapter allows users to set the microscope's eyepiece height to achieve a good ergonomic neutral posture. Adjustability eliminates raised shoulders and bent necks, which can cause discomfort.</p>




Laboratory Ergonomics	Features/Options	Image	Applications
<p>Extended Eyetube Binocular</p> <p>Source:  <a href="http://www.bayoptical.com/extendedeyetube.html">http://www.bayoptical.com/extendedeyetube.html</a></p>	<ul style="list-style-type: none"> <li>- 90MM extension</li> <li>- Available as a new item or as a retrofit of existing binocular head (depending on the binocular model)</li> <li>- Can be combined with the ErgoAdapter and an Optical Wedge for greatly enhanced workstation flexibility</li> </ul>		<p>Extended Eyetube Binocular is an ergonomic device for use with quality stereomicroscopes. Its longer-than-standard eyetube length enables the microscope user to comfortably move back from the workstation, helping to relieve postural fatigue. Ideal for most microscopes, especially those used with machines, bonders, heated stages, solder stations, under hoods or anywhere operators are otherwise forced into a cramped viewing posture.</p>
<p>Optical Wedge</p> <p>Source:  <a href="http://www.bayoptical.com/opticalwedge.html">http://www.bayoptical.com/opticalwedge.html</a></p>	<ul style="list-style-type: none"> <li>- Easily installed due to self-aligning flanges</li> <li>- Offers a 0° or 90° in-line viewing angle</li> <li>- Maintains parfocality and optical alignments at any inclination angle</li> </ul>		<p>Adds new flexibility to microscope workstations. 0°/90° Optical Wedge is positioned between the binocular head and the optics carrier. It complements existing 45° binoculars, providing two additional viewing angles. Costs less than half of a straight tube binocular and offers three times the flexibility.</p>
<p>ISIS Ergonomic Eyepiece</p> <p>Source:  <a href="http://www.optimaxonline.com/newsdetails.php?newsId=23">http://www.optimaxonline.com/newsdetails.php?newsId=23</a></p>	<ul style="list-style-type: none"> <li>- Fits most brands of conventional binocular microscopes</li> <li>- Expanded-Pupil technology expands field of vision</li> <li>- Enables wearing of prescription spectacles</li> <li>- Eliminates direct contact with the eyepieces</li> <li>- Removes optical the disturbances known as mouches volantes (floaters)</li> <li>- Far less costly than replacing entire instrument</li> </ul>		<p>Vision Engineering's ISIS microscope eyepieces are patented 'Expanded-Pupil' eyepieces designed to replace the conventional eyepieces of well known, other brand instruments. ISIS increases radial and axial head freedom by at least 8 times that of conventional eyepieces and greatly improves the visibility of optical information across the entire field. A unique eye distance of 38mm allows glasses to be worn in complete comfort, with common weaknesses (long and short sight astigmatism) no longer affecting user efficiency. Designed to fit most microscopes with standard diameter eyepiece tubes (23mm or 30mm) and can be used with many microscopes, including Kyowa, Leica, Meiji, Motic, Nikon, Olympus, Wild, Zeiss.</p>

Laboratory Ergonomics	Features/Options	Image	Applications
<p>Fluorescent Ring Light System</p> <p>Source:  <a href="http://www.microscopeworld-professional.com/detail.aspx?ID=251">http://www.microscopeworld-professional.com/detail.aspx?ID=251</a></p>	<ul style="list-style-type: none"> <li>- Your choice of 110v or 220v power supply</li> <li>- Inside diameter of ring is 60mm</li> <li>- Lamp life 1,500 hours</li> <li>- Color temperature 6500K</li> <li>-</li> </ul>		<p>Lightweight, compact, and rugged, this ring illuminator provides 360 degrees of cool, white, shadow-free illumination. For use with EM Series stereo microscopes (adapter required).</p>
<p>Microscope Video Systems</p> <p>Sources:</p> <ul style="list-style-type: none"> <li>- <a href="http://www.scientificdevice.com">www.scientificdevice.com</a></li> <li>- <a href="http://www.usbdigitalmicroscope.com">www.usbdigitalmicroscope.com</a></li> <li>- <a href="http://www.office-products.cn/trade/nanonist-2-0-mega-digital-microscope-eyepiece-11375">www.office-products.cn/trade/nanonist-2-0-mega-digital-microscope-eyepiece-11375</a></li> <li>- <a href="http://www.microscopesfromnightengale.com">www.microscopesfromnightengale.com</a></li> </ul>	<ul style="list-style-type: none"> <li>- Easy to attach to a microscope</li> <li>- Variety of magnifications available</li> <li>- Provides wider viewing area and increased working distance</li> </ul>		<p>Camera systems provide remote viewing of microscope image on computer monitor to eliminate awkward head and neck postures associated with viewing through an eyepiece.</p>
<p>Meiji Graduated Mechanical Stage MA564</p> <p>Source:  <a href="http://www.microscopeworld-professional.com/detail.aspx?ID=65">http://www.microscopeworld-professional.com/detail.aspx?ID=65</a></p>			<p>Graduated mechanical stage for Meiji microscope to reduce awkward reaches required to position lab tissues and cells for observation.</p>
<p>Microdesk</p> <p>Source:  <a href="http://www.backdesigns.com/store/Regular-Microdesk-P555C66.aspx">http://www.backdesigns.com/store/Regular-Microdesk-P555C66.aspx</a></p>	<ul style="list-style-type: none"> <li>- Transparent acrylic</li> <li>- Adjustable height and slope</li> <li>- Two sizes available</li> <li>- Removable bottom lip accommodates large documents and media (e.g., blueprints, x-rays).</li> <li>- Easily lifted and stored when not in use.</li> </ul>		<p>The Microdesk is an adjustable sloped writing platform for the keyboard work position. It enables easy access to the keyboard while working with papers and books. Prevents reaching and twisting as your papers are placed directly in front of you. Can be used to "stack" work to keep materials within close reach zone.</p>



Laboratory Ergonomics	Features/Options	Image	Applications
Pipettes and Accessories			
<p>Rainin Pipet Lite LTS</p>	<ul style="list-style-type: none"> <li>- Reduces pipetting force by 70%</li> <li>- Magnetic assist locates and holds zero position</li> <li>- LTS tip ejection system reduces tip ejection force</li> </ul>		<p>Pipet Lite reduces total force by 70% .Magnetic assist helps locate and hold zero position for improved results. LTS® LiteTouch Tip Ejection for ultralight tip ejection. Reduces your risk of RSIs. Uses LTS® tips only</p>
<p>Rainin EDP3 Plus Advanced Electronic Pipette (Single and Multichannels available)</p> <p>Source:  <a href="http://www.rainin.com/products/product_list.asp?class=46">http://www.rainin.com/products/product_list.asp?class=46</a></p>	<ul style="list-style-type: none"> <li>- Reduces pipetting force by 95%</li> <li>- Seven pipetting cycles</li> <li>- LTS tip ejection system reduces tip ejection force</li> </ul>		<p>Reduce pipetting forces by 95%. Advanced Ergonomic Pipette  Pipet, multidispense, titrate, dilute, volume sequence, GLP, and seven options.LTS® tips for ultralight tip ejection. Includes Li-Ion battery, 3000 cycles per charge. Order charging devices separately</p>
<p>Vistalab Ovation Manual and Electronic Single and Multiple Channel Pipettes</p> <p>Source:  <a href="http://www.vistalab.com/OvationProductsMenu.asp">http://www.vistalab.com/OvationProductsMenu.asp</a></p>	<ul style="list-style-type: none"> <li>- Low-force, form-fitting plunger manually controls aspiration and dispensing speeds with minimal effort</li> <li>- Light tip ejection force</li> <li>- Adjustable finger hook</li> <li>- Disposable battery</li> </ul>		<p>Unique design of Vistalab Ovation reduces shoulder elevation, wrist deviation and wrist flexion and extension.</p>
<p>Pipette Carousels</p> <p>Source:  <a href="http://www.biohit.com/view/products.asp?document_id=1029&amp;cat_id=276">http://www.biohit.com/view/products.asp?document_id=1029&amp;cat_id=276</a>  <a href="http://www.biohit.com/view/products.asp?document_id=1029&amp;cat_id=276">http://www.biohit.com/view/products.asp?document_id=1029&amp;cat_id=276</a></p>	<ul style="list-style-type: none"> <li>- Carousel for 6 pipettes in a minimal space for storing all Research pipettes (single and multi-channel)</li> <li>- Carousels available for variety of models of pipettes</li> </ul>		<p>The rotating pipette stand makes for an organized workplace. The new kind of attachment means that the pipette is no longer held at the bottom, ruling out potential contamination of a pipette by it being put on a stand.</p>




Laboratory Ergonomics	Features/Options	Image	Applications
<p>FinnPipette Novus Electronic Multi Channel Pipette</p> <p>Source:  <a href="http://www.thermo.com/com/cda/product/detail/1,1055,10120040,00.html">http://www.thermo.com/com/cda/product/detail/1,1055,10120040,00.html</a></p>	<ul style="list-style-type: none"> <li>- Backlit display</li> <li>- Trigger action buttons allow dispensing with index fingers reducing thumb motion</li> <li>- Adjustable finger rest</li> <li>- Lightweight</li> <li>- 10 pipetting functions</li> <li>- 9 speeds</li> <li>- 6 language interface</li> <li>-</li> </ul>		<p>Finnpipette Novus is the first electronic pipette to have backlit technology. It is lightweight and has an adjustable finger rest for use with either hand. Its multiple pipetting functions and multiple speeds make it one of the most versatile electronic pipettes on the market.</p>
<p>Magnifier Lamps # 45-3357</p> <p>Source:  <a href="http://www.sciencelab.com/page/S/PVAR/10175/45-3357">http://www.sciencelab.com/page/S/PVAR/10175/45-3357</a></p>	<ul style="list-style-type: none"> <li>- 2X magnification, 36 in. reach</li> <li>- Convenient front handle</li> <li>- Spring counterbalance</li> <li>- Control knobs lock magnifier in place</li> <li>- Large viewing area, 4 x 21/2 in.</li> <li>- Incandescent 60w bulb</li> <li>- 120V 60Hz, 0.5A</li> <li>- UL listed and CSA approved</li> </ul>		<p>Adjustable position light magnifier can reduce awkward postures and eye strain when viewing laboratory tasks that require magnification.</p>
<p>Laboratory Scissor Jack - Stainless Steel, Plate Size 6"x6" - 10-9003</p> <p>Laboratory Scissor Jack - Stainless Steel, Plate Size 8"x8" - 10-9004</p> <p>Source:  <a href="http://www.sciencelab.com/page/S/PVAR/10-90034">http://www.sciencelab.com/page/S/PVAR/10-90034</a></p>	<ul style="list-style-type: none"> <li>- Laboratory Jack 10-9003 collapsed height is 3" and expanded height is 10".</li> <li>- Laboratory jack 10-9004 collapsed height is 3.5" and expanded height is 10.25".</li> </ul>		<p>Laboratory scissor jacks constructed of high quality stainless steel to raise equipment height. Supports up to 200 lbs. Two plate sizes.</p>
<p>Test Tube Clamps</p> <p>Source:  <a href="http://www.indigo.com/science-supplies/gph-science-supply/lab-support-stand.html">http://www.indigo.com/science-supplies/gph-science-supply/lab-support-stand.html</a></p>	<ul style="list-style-type: none"> <li>- Clamps can support test tubes and flasks</li> <li>- Support rings hold beakers &amp; flasks</li> <li>-</li> </ul>		<p>Clamps with stands that hold and position test tubes or cylinders to minimize static postures associated with holding them for extended periods of time</p>






Laboratory Ergonomics	Features/Options	Image	Applications
<p>Flexible Gooseneck Arms with Electrode Holder</p> <p>Source:  <a href="http://www.belart.com/shop/183152322-flexible-with-weighted-base-p-183152322.html?cPath=19">http://www.belart.com/shop/183152322-flexible-with-weighted-base-p-183152322.html?cPath=19</a></p>	<ul style="list-style-type: none"> <li>- 46cm (18") long x 25mm (1") in diameter arm</li> <li>- One end has a 12.7mm (1/2") diameter stainless steel pin for attachment to a clamp or base.</li> <li>- Attach electrode holder to other end</li> <li>- Molded of chemical and corrosion resistant acetal copolymer links with ball-and-socket joints</li> <li>- Base and electrode holder sold separately• Chemical resistant</li> </ul>		<p>Flexible gooseneck arm and base with electrode holder can be used to position a Ph meter or similar tool and eliminate prolonged gripping of tools.</p>
<b>Bottles and Dispensers</b>			
<p>Plastic Wash Bottles (Squeeze Bottles)</p> <p>Source:  <a href="http://www.daigger.com/catalog/product?deptId=&amp;prodId=2757D">http://www.daigger.com/catalog/product?deptId=&amp;prodId=2757D</a></p>	<ul style="list-style-type: none"> <li>- Low density polyethylene with leakproof screw cap</li> <li>- 1-pc. cap/spout</li> </ul>		<p>Transfer small amounts of unmeasured liquids in place of test tubes or beakers.</p>
<p>Graduated Cylinder with Handle</p> <p>Source:  <a href="http://www.usplastic.com/catalog/product.asp?catalog%5Fname=USPlastic&amp;category%5Fname=14&amp;product%5Fid=12826&amp;cookie%5Ftest=1">http://www.usplastic.com/catalog/product.asp?catalog%5Fname=USPlastic&amp;category%5Fname=14&amp;product%5Fid=12826&amp;cookie%5Ftest=1</a></p>	<ul style="list-style-type: none"> <li>- Standard polypropylene</li> <li>- Double scale and octagonal base</li> <li>- Steam autoclavable at 250°F</li> </ul>		<p>The generously sized handle makes pouring from a tall cylinder safer and more convenient</p>
<p>3181 Bottle Jacket</p> <p>Source:  <a href="http://www.labsafety.com/search/bottle+jacket/35585/3181/?isredirect=true">http://www.labsafety.com/search/bottle+jacket/35585/3181/?isredirect=true</a></p>	<ul style="list-style-type: none"> <li>- Tough polyethylene shell encapsulates a 4L or 1-gal. bottle</li> <li>- Two-piece jacket with unique threading design firmly connects top and bottom</li> <li>- Blocks out light to reduce solvent degradation.</li> </ul>		<p>Easy grip handle improves ability to grip and pour gallon or 4-liter bottles of liquids. Reduces spill hazards in case of breakage.</p>



Laboratory Ergonomics	Features/Options	Image	Applications
<p>Brinkmann Bottletop Dispensers</p> <p>Source:  <a href="http://www.daigger.com/catalog/product?deptId=&amp;prodId=20565&amp;q=bottle+dispensers">http://www.daigger.com/catalog/product?deptId=&amp;prodId=20565&amp;q=bottle+dispensers</a></p>	<ul style="list-style-type: none"> <li>- Precise volume settings</li> <li>- Filling tube adjusts to fit reagent bottle</li> <li>- Dispensers rotate 360°</li> <li>- Calibration certificate</li> </ul>		<p>The Brinkmann Bottletop Dispensers make reagent dispensing as quick and easy as possible. Dispensers are available in fixed-volume or adjustable volume sizes to cover the range of 0.5 to 100 mL.</p>
<p>Two Handle Dispensing Jugs</p> <p>1026 Nalgene Safety Dispensing jug</p> <p>Source:  <a href="http://www.labsafety.com/search/jug/24529963/">http://www.labsafety.com/search/jug/24529963/</a></p>	<ul style="list-style-type: none"> <li>- Low-density polyethylene jug includes top and bottom safety handle</li> <li>- Guaranteed leakproof</li> <li>- Permanently attached pour spout reduces splash</li> <li>- 83mm polypropylene closure</li> <li>- Not autoclavable.</li> </ul>		<p>Double-handle design makes carrying and pouring easier and safer to accomplish.</p>
<p>Bottom Dispensing Carboy</p> <p><a href="http://www.daigger.com">www.daigger.com</a>  <a href="http://www.usplastic.com">www.usplastic.com</a></p>	<ul style="list-style-type: none"> <li>- Bottom dispenser</li> <li>- Rounded octagonal shape increases strength</li> <li>- Offset neck makes pouring easier.</li> </ul>		<p>Bottom dispensing carboy used to dispense large volumes of liquids from bottom of container. Eliminates need to lift bottle and pour or twist off cap or lid.</p>



Laboratory Ergonomics	Features/Options	Image	Applications
<p><b>Chairs, Footrests and Mats</b></p> <p>Humanscale Freedom Saddle Stool</p> <p>Source:  <a href="http://www.humanscale.com/products/saddle_seat.cfm">http://www.humanscale.com/products/saddle_seat.cfm</a></p>	<ul style="list-style-type: none"> <li>- Uniquely contoured triangular cushion allows multiple sitting positions for ultimate comfort</li> <li>- Height adjustable pneumatic cylinder accommodates a wide range of users and tasks.</li> <li>- Optional foot ring and high cylinders for specific applications.</li> <li>- Available in a wide variety of textiles, including wipeable vinyls for lab, cleanroom and healthcare use.</li> </ul>	 <p style="text-align: center;"><i>freedom</i> Saddle Seat</p> <p style="text-align: center; font-size: small;">click here to view Saddle Seat without footring option</p>	<p>The Saddle seat, part of Humanscale's award-winning Freedom seating line, is the most comfortable and versatile stool ever made. Its unique, triangular cushion encourages users to sit in a "saddle" posture, which lowers the thighs, opens up the hips and puts the spine into a healthy lordotic curve. The "saddle" posture also reduces pressure points for long-term comfort and allows users to get closer to their work. It's simply one of the most comfortable and healthy ways to sit!</p>
<p>Bambach Stool</p> <p>Source:  <a href="http://www.backdesigns.com/Bambach-Saddle-Chair-P566C158.aspx">http://www.backdesigns.com/Bambach-Saddle-Chair-P566C158.aspx</a></p>	<ul style="list-style-type: none"> <li>- Seat tilts 5° backward and 10° forward for comfort.</li> <li>- 27" base is provided with tall lifts used with a footring for greater stability at no additional charge.</li> <li>- Oversize 60mm hooded casters for smooth rolling on carpet or hard floors.</li> <li>- Optional braking, locking, and industrial casters also available.</li> </ul>		<p>The Bambach Saddle Seat is shaped like a saddle to enable the body to maintain ideal posture when sitting. There is no need for direct lumbar support because the correct lumbar curve is maintained naturally rather than artificially.</p>

Laboratory Ergonomics	Features/Options	Image	Applications
<p>Hag Capisco Laboratory chair</p> <p>Source:  <a href="http://www.ergodepot.com/HAG_Ultrasound_Chair_p/8106u.htm">http://www.ergodepot.com/HAG_Ultrasound_Chair_p/8106u.htm</a></p>	<ul style="list-style-type: none"> <li>- Non-traditional saddle seat design. Flat seat is available at no additional charge.</li> <li>- Adjustments: seat depth, seat height, back height, tilt tension backward, and three-position tilt lock.</li> <li>- Provides full trunk and arm support when seated backwards.</li> <li>- Aluminum base with curved/arched footring</li> </ul>		<p>HAG's award-winning Capisco chair can be modified for ultrasound, radiology and laboratory applications. The tall 265mm lift offers nearly ten inches of seat height variation. And the footring (or optional Step-Up) allows for infinite positions for your feet. Ideal for sit/stand applications, the chair encourages creative, alternate sitting positions - forwards, backwards or sideways.</p>
<p>Labster BIMOS Chair</p> <p>Source:  <a href="http://www.industrysearch.com.au/News/New_Bimos_Labster_Worlds_First_Pure_Lab_Chair-35827">http://www.industrysearch.com.au/News/New_Bimos_Labster_Worlds_First_Pure_Lab_Chair-35827</a></p>	<ul style="list-style-type: none"> <li>- Clean design with no sharp edges or crevices</li> <li>- Covered star base for stability and easy cleaning</li> <li>- Height adjustable footrest</li> <li>- Compact footprint</li> <li>- Seat tilt</li> <li>- Tilting molded backrest</li> <li>-</li> </ul>		<p>The <b>Labster</b> from BIMOS Germany is specifically designed for <b>laboratory work</b>, with no sharp edges or crevices and a closed base for easy cleaning. Compact design provides small footprint for small spaces.</p> <p>Labster's special "work function" backrest with specially moulded polyurethane features a unique ventral support. The seat unit <b>flexibly</b> adapts to the corresponding position of the body and permits ergonomic work in postures that are typical for labs, where operators often spend long hours with the upper part of their body bent forwards, for example when pipetting or working at a microscope.</p>
<p>Adjustable footrest SV-973-700</p> <p>Source:  <a href="http://www.scientificvisions.com/catalog/product_info.php?cPath=26&amp;products_id=33">http://www.scientificvisions.com/catalog/product_info.php?cPath=26&amp;products_id=33</a></p>	<ul style="list-style-type: none"> <li>- Adjustable</li> <li>- Slide Resistant Rubber Feet</li> <li>- Easily Cleaned</li> <li>- Non-skid Foot Surface</li> <li>-</li> </ul>		<p>Adjustable, tall, heavy duty footrest provides support for feet when sitting at lab bench or using safety cabinet.</p>

Laboratory Ergonomics	Features/Options	Image	Applications
<p>ISE Steel Frame Footrest</p> <p>Source:  <a href="http://airtechcorp.stores.yahoo.net/qfstlisestfr.html">http://airtechcorp.stores.yahoo.net/qfstlisestfr.html</a></p>	<ul style="list-style-type: none"> <li>- 23.5" w x 13.5"d x 17" h</li> <li>- Adjusts in 1" increments between 4-14" in height</li> <li>- 3 tilt angle adjustments</li> <li>- Steel platform on non-skid pads</li> </ul>		<p>Very sturdy footrest can be used to support feet when working at lab bench or in safety cabinet or hood.</p>
<p>BioFit Industrial Footrest</p> <p>Source:  <a href="http://www.labsafety.com/BIO-FIT-Industrial-Footrests_24548808/">http://www.labsafety.com/BIO-FIT-Industrial-Footrests_24548808/</a></p>	<ul style="list-style-type: none"> <li>- Adjust height in 1 in. increments from 3-11 in.</li> <li>- Fixed 10° slope</li> <li>- Nonslip rubber matting coating</li> <li>- Durable steel tubing frame with powder coat finish</li> </ul>		<p>Industrial grade footrest that can support feet when working at a lab bench or safety cabinet.</p>
	<ul style="list-style-type: none"> <li>- Height adjustment in 1" increments from 6" to 15" above the floor</li> <li>- Fluid rocking motion automatically adjusts with changes in body position.</li> <li>- Heavy gauge steel frame</li> <li>- Large 23"w x 11.25"d platform</li> </ul>		<p>Adjustable, tall, and heavy duty footrest provides support for feet when sitting at lab bench or using safety cabinet.</p>
<p>Ortho Mat Floor Mats</p> <p>Source:  <a href="http://airtechcorp.stores.yahoo.net/protormat.html">http://airtechcorp.stores.yahoo.net/protormat.html</a></p>	<ul style="list-style-type: none"> <li>- Closed cell mat eliminates absorption of liquids, fluids and airborne pathogens</li> <li>- Infection control "antimicrobial" prevents mold, mildew and bacterial growth</li> <li>- Light weight, easy to clean. disinfectants.</li> <li>- Custom sizes available</li> </ul>		<p>Floor mats designed to reduce back, leg, foot, and ankle fatigue when standing for prolonged periods to complete lab tasks.</p>
<p>Dawg Vinyl Floor Mat 488</p> <p>Source:  <a href="http://www.dawginc.com/floor-safety/anti-fatigue/vinyl-mat.php">http://www.dawginc.com/floor-safety/anti-fatigue/vinyl-mat.php</a></p>	<ul style="list-style-type: none"> <li>- Vinyl non-porous surface with Nitricell sponge base</li> <li>- Anti-microbial</li> <li>- Abrasion and chemical resistance</li> <li>-</li> </ul>		<p>Floor mats designed to reduce back, leg, foot, and ankle fatigue when standing for prolonged periods to complete lab tasks.</p>

Laboratory Ergonomics	Features/Options	Image	Applications
<p>SATech SmartCell Anti-Fatigue Mats and Flooring</p> <p>Source:  <a href="http://www.satechinc.com/products_mats.htm">http://www.satechinc.com/products_mats.htm</a></p>	<ul style="list-style-type: none"> <li>- Durable natural rubber</li> <li>- Matrix of hollow cylinders that function like spring</li> <li>-</li> </ul>		<p>SmartCells® anti-fatigue mats enable workers to stand comfortably for long periods while minimizing the fatigue and pain associated with standing work.</p>
<p>Countertop Extender</p> <p>Source:  <a href="http://cndproducts.com/labcare.php">http://cndproducts.com/labcare.php</a></p>	<ul style="list-style-type: none"> <li>- Dual slide out work surfaces</li> <li>- Shelves then store neatly under the raised work surface when not in use</li> <li>- Non-skid, anti-vibration feet</li> <li>- Can be cleaned with disinfectant</li> </ul>		<p>Countertop Extender instantly adds valuable work space to any counter. It has dual slide out work surfaces to allow extra work space on any bench. The two spring-assisted shelves then store neatly under the raised work surface when not in use.</p> <p>The rugged construction can support even heavy bench top equipment. The non-skid, anti-vibration feet will keep the countertop extender in place.</p>